

Hip Hold Variable Angle – *Directions for Use*

1) Place the patient on the operating table with the affected hip upwards.

2) *Posterior*

Mount the vertical posterior support on the operating table.
Adjust the support bars so the pad is against the patient's sacrum.
Lock the vertical support bar to the operating table.

3) *Anterior*

Mount the vertical anterior support on the operating table.
Adjust both anterior supports so they are positioned against the patient's iliac crest.

Lock the vertical support bar to the operating table.

Push the lower horizontal support lightly against the operating table.
Tighten the Bakelite knob.

Adjust the patient so the pelvis is perpendicular to the table.
Lock the upper horizontal support with the locking screw.

4) *Posterior*

Press the posterior horizontal support against the patient to get stable fixation.
Lock the Bakelite knob.

Notes

- If there is a risk of skin necrosis, additional cushioning should be applied to the support pad. The gel pad is ideal.
- Check that the support cushion and lock handle are fixed in the support arm.
- If any part of the Hip Fix shows signs of wear, it is essential that the part be replaced so that optimal results and absolute stability can be maintained.
- Cease use immediately if any parts appear faulty and contact the manufacturer.
- The Bakelite Knob should be tightened by HAND using moderate force.
- Check for any wear of the operating table attachments so that they do not negatively affect the use and mounting of the Hip Hold.
- Weight limit of patient 150kg.




fig 1



fig 2




fig 3




fig 4